

PRX-PLUS®

THE EVOLUTION OF NEEDLE-FREE BIOREVITALIZATION



Dear Partner,

We are pleased to introduce **WiQo PRX-PLUS**, the next-generation professional skin booster designed to meet the growing demand for visible, natural, and immediate results – without invasiveness or downtime.

Building on the innovation of the renowned PRX-T33, PRX-PLUS enhances its benefits and is particularly suited for patients seeking **comfort, safety, and practicality**.

Key Features of PRX-PLUS

- **Patented Formula:** Combines Ammonium Trichloroacetate (ATCA) – a gentler form of TCA – with exclusive Homocysteic Acid (HoA), plus citric and phytic acids.
- **Mechanism of Action:** Stimulates deep dermal regeneration, improving skin radiance, firmness, texture, and tone – with no visible peeling.
- **Fast & Comfortable Protocol:** Painless, needle-free, and no downtime.
- **Safe & Versatile:** Usable all year round, including on dark phototypes and in summer, with no significant side effects.
- **Longer Shelf-Life:** 30 months, compared to 15 for PRX-T33.
- **Easier Storage:** No refrigeration required, simplifying logistics and management.

Clinically Proven Results

After treatment:

- Wrinkle depth reduced by up to **40%**
- Decrease in **microroughness and enlarged pores**
- **Sebum regulation**
- Enhanced **evenness and brightness of complexion**

What Sets PRX-PLUS Apart

Although part of the same needle-free family as PRX-T33, PRX-PLUS is an evolution, thanks to:

- **Greater Tolerability:** ATCA is less aggressive than 33% TCA, providing a more comfortable experience.
- **Exclusive Homocysteic Acid (HoA):** Enhances regenerative effect without increasing irritation.
- **Longer Shelf Life & No Cold Chain:** 30 months, no refrigeration.
- **Sensitive Skin Friendly:** Great for patients with reactive skin or treatment phobia.
- **Fully Seasonal-Compatible:** Truly no downtime, even in summer.
- **Up to 20 Treatments per Box**

PRX-T33	PRX-PLUS
Deep photoaging	Dull, lackluster skin
Melasma and dark spots	Uneven texture, enlarged pores
Acne scars	Superficial roughness
Advanced wrinkles and laxity	Early signs of aging, loss of tone
Patients willing to undergo more intense treatment	Ideal for patients who fear needles, pain, or downtime
	Ideal for maintenance and anti-aging prevention
	Suitable in contexts requiring logistical simplicity and long shelf-life

Ideal Patients

PRX-T33: Motivated, correction-oriented patients needing deeper stimulation.

PRX-PLUS: Younger or mildly affected patients seeking a gentle but effective approach – perfect for summer and darker skin types.

Conclusion

WIQO PRX-PLUS expands the range of needle-free treatments, complementing and enhancing PRX-T33. It is not a replacement, but a **complementary option**, allowing practitioners to **tailor the intensity of dermal stimulation** based on the clinical picture and patient expectations.

Integrating both protocols into your practice enables a personalized approach for each patient, offering **scientifically validated, safe, and high-performance solutions**—with PRX-PLUS providing the added benefit of simpler and more convenient product management, especially from a logistical standpoint.

/Before&After

before



after 3 sessions



Skin Laxity: The tightening effect visibly improves nasolabial folds, marionette lines, and overall facial contours.

before


Wrinkles: Smoothing effect reduces wrinkle depth by up to 40%.



after 4 sessions



before



after 5 sessions



Uneven Skin Tone: The brightening effect promotes a more uniform pigmentation.

/Treatment Protocol

PRX-PLUS TREATMENT: THE NEEDLE FREE SKIN BOOSTER IN 3 SIMPLE STEPS



1\ Prepare



2\ Stimulate



3\ Nourish

Full cycle: 4 sessions on average

Session duration: 20–25 min (face) / 30 min (face + neck)

Product amount: 1 ml (face) / 1.5 ml (face + neck)

Session frequency: Every 14 days

Maintenance sessions: Every 45 days (post-cycle)

Second cycle: After 6 months